



# LIGNIN

## pH Electrode

### QUICK PREP

#### Fill with $\text{KNO}_3$

Slide cover-sleeve down to expose fill-hole  
Fill chamber with 10%  $\text{KNO}_3$  or 3M  $\text{KCl}^*$   
solution

Press and release cap to allow junction to fill  
with solution

Top off chamber if necessary

#### Soak overnight in pH 4.01 buffer solution

Drain and refill with  $\text{KNO}_3$  or  $\text{KCl}^*$

Rinse

### OPERATION

#### Calibration

Perform a two-point calibration with fresh  
buffer solutions

#### Use

Rinse with de-ionized or distilled water  
between measurements

Press cap to flush reference junction with  
filling solution after problem samples (and  
refill reservoir if empty)

#### Storage

Soak in pH electrode storage solution or 4.01  
buffer solution

\*3M  $\text{KCl}$  solution = 225 g  $\text{KCl}$  per 1 L  $\text{H}_2\text{O}$



## MAINTENANCE

### Disassembly

- Unscrew cap and remove black retaining ring
- Unscrew cap-body and slide it (along with the spring) back from the main assembly
- Push the internal element out of the body to expose the tapered reference junction

### Cleaning

- Wash the junction with dish-soap and warm tap water (do not use solvents)
- Rinse with de-ionized or distilled water

### Assembly

- Gently pull cable to re-seat the junction inside the electrode body
- Replace the spring and cap-body and re-thread the cap-body for a snug but not over-tightened fit
- Replace the black retaining ring and route the cable through its groove
- Screw on the cap

### Long-Term Storage

- Drain reservoir
- Flush with de-ionized or distilled water
- Ensure reservoir and junction are dry and cover electrode with protective guard